


Andy French, Tony Ayres, Louisa Lintern Jones. Department of Physics, Winchester College. 19/20 December 2016



**Isaac Newton** proposed a **Law of Universal Gravitation** to calculate the *attractive* force between point masses m and M, separated by a distance r

$$F = \frac{GMm}{r^2}$$

Isaac Newton 1642-1727


The constant G is:

$$G = 6.67384 \times 10^{-11} \,\mathrm{m^3 kg^{-1} s^{-2}}$$




Henry Cavendish 1731-1810

To measure *G* experimentally on Earth is *difficult*, since gravitational forces only become important on a human scale when the masses are *planet sized*.



Henry Cavendish devised an ingenious experiment to overcome this problem. His equipment was based upon a **torsional pendulum**. Amazingly, the tiny forces resulting from the gravitational attraction of lead spheres could result in a *measurable* effect upon the oscillation of the pendulum.



We can now combine the

expressions to find G

Measuring G via the Cavendish experiment

Balance the torsion force (the twist) on the wire with the torque resulting from the gravitational attraction of masses M and m

The torsion constant  $\tau$  can be found by measuring the period T of small oscillations of the pendulum.

$$I = 2 \times m \left(\frac{1}{2}L\right)^2 = \frac{1}{2}mL^2$$
$$\frac{T^2}{4\pi^2} = \frac{\frac{1}{2}mL^2}{\tau}$$
$$= \frac{2\pi^2 mL^2}{\tau}$$

 $T = 2\pi \sqrt{\frac{I}{\tau}} \leftarrow$ 

 $- T^{2}$ 

For a mass on a spring of stiffness *k* 

 $T = 2\pi \sqrt{\frac{m}{k}}$ 

Angle in

radians

 $\tau\theta = L \times \frac{GmM}{r^2}$ 

 $\therefore G = \frac{\tau r^2 \theta}{LmM}$ 

Moment of inertia I of the pendulum about the wire axis

 $G = \frac{2\pi^2 L r^2 \theta}{2\pi^2 L r^2 \theta}$ 


 $G = \frac{2\pi^2 m L^2}{T^2} \times \frac{r^2 \theta}{LmM}$ 

In the original experiment m = 0.73kg, M = 158kg L = 1.8m, r = (230 - 4.1)mm P = 875.3s,  $\frac{1}{2}L\theta \approx 4.1$ mm


Cavendish measured

$$G = 6.74 \times 10^{-11} \text{ m}^3 \text{kg}^{-1} \text{s}^{-2}$$
$$[G_{\text{modern}} = 6.67 \times 10^{-11} \text{ m}^3 \text{kg}^{-1} \text{s}^{-2}]$$

At Winchester, we have a Cavendish Experiment kit supplied by LD Didactic GmbH




Small lead masses, fixed to torsional pendulum m = 0.015 kg



position – as in our case given we have t platform installed in Laboratory P4.




The laser is reflected off a mirror, which rotates with the pendulum. This causes a vertical stripe to slowly track back and forth across the opposing wall of the laboratory. The dynamics of these oscillations are used to determine G.



Once the pendulum was aligned, measurements were taken, every 30s, from the *centre* of the stripe along the tape measure scale.

A measuring tape was affixed to the wall such that the stripe tracks along it.





These screws are used to lower metal 'baskets' which lock the pendulum in place. The baskets *must* be raised when the equipment is transported.




Horizontal alignment screws

Once the lowered baskets have released the pendulum, it is often discovered that the pendulum will 'crash' into the glass windows. i.e. it has too much angular momentum and the equilibrium position is not central. The top thumb – screw enables the wire and mirror to be manually rotated to correct for this.


This process of alignment can be time consuming and turns of the screw should be subtle! This screw seems to relate to how the wire is fixed in place. Turning it has broken the kit in the past....





Diagrams from LD Didactic Instructions P1.1.3.1

is rotated till it just touches the glass panel



The offset angle of the laser will mean the equilibrium positions are offset from N

We actually measured to a spot just above the 'window' rather than the mirror, as the range-finder tended to pick up the back wall otherwise.

i.e. add 29/2 mm to the range measurements.




A **laser range-finder** was used to determine the distances between the equipment and the wall.

$$L_0 = 6820 \text{mm}$$

The distance from the perpendicular point N to the **equilibrium of Position A** was measured to be

$$L_1 = 1292 \,\mathrm{mm}$$

**Errors** of about 10mm (0.8%) due to wobble of the laser when button pressed, misalignment of walls etc




Once the system was aligned, and it was confirmed that oscillations were contained within the pendulum housing, (i.e. no crashing against the interior windows) laser stripe positions along the tape were recorded every 30s.

After four oscillations, Position A was changed to Position B, and another three complete oscillations were recorded.

This amounted to an experimental run time of about 60 minutes.


Data was recorded directly into **Microsoft Excel** and plotted as the data was collected.



Note the equilibrium stripe location for position A will be **inferred** via curve fitting (rather than wait many hours!)

This could be a source of error in both **period** and **AB equilibrium offset** 






After several hours following the cessation of measurements recorded in **Position B**, the following measurement was photographed. This is assumed to be the **Position B equilibrium reading**:

## 4032mm


The data was then imported into **MATLAB** for further analysis.

- **1. Cubic spline interpolation** to smooth between the data points.
- 2. Peak (and trough) finding routine to determine the peak and trough times and positions for each A and B oscillation.
- 3. The difference between peak times were also used to calculate the **period** of torsional oscillations, which was calculated to be **616s.**
- 4. Curve fitting to estimate the longer term oscillation, and predict the final equilibrium positions.
- 5. Calculation of *G* (see later slide for details)
- 6. Automatic graph plot and .png file creation.





## Curve fitting to determine the equilibrium position



$$\frac{x_0 - x_1}{x_2 - x_3} = \frac{A + Ae^{-\frac{1}{2}\gamma T}}{Ae^{-\gamma T} + Ae^{-\frac{3}{2}\gamma T}}$$
$$\frac{x_0 - x_1}{x_2 - x_3} = \frac{1 + e^{-\frac{1}{2}\gamma T}}{e^{-\gamma T} \left(1 + e^{-\frac{1}{2}\gamma T}\right)} = e^{\gamma T}$$
$$\therefore \gamma = \frac{1}{T} \ln\left(\frac{x_0 - x_1}{x_2 - x_3}\right)$$

$$x(t) = Ae^{-\gamma(t-t_{\max})} \cos\left(\frac{2\pi \left(t-t_{\max}\right)}{T}\right) + x_{\infty}$$
$$x_{0} = A + x_{\infty}$$
$$x_{1} = -Ae^{-\frac{1}{2}\gamma T} + x_{\infty}$$
$$x_{2} = Ae^{-\gamma T} + x_{\infty}$$
$$x_{3} = -Ae^{-\frac{3}{2}\gamma T} + x_{\infty}$$

$$\frac{x_0 - x_1}{x_2 - x_3} = \frac{A + Ae^{-\frac{1}{2}\gamma T}}{Ae^{-\gamma T} + Ae^{-\frac{3}{2}\gamma T}}$$
$$x_0 - x_1 = A\left(1 + e^{-\frac{1}{2}\gamma T}\right)$$
$$\therefore A = \frac{x_0 - x_1}{1 + e^{-\frac{1}{2}\gamma T}}$$
$$\therefore x_\infty = x_0 - A$$





Diagram and parameters from LD Didactic literature

We measured this M = 1.4966 kg R = 32 mm m = 0.015 kg r = 6.9 mmLarge sphere Small sphere

d = 50mm b = 46.5mm

> The French language version (supplied with the kit!) states 46.5mm whereas the English language version states 50mm....



Antitorque moment of the "second" lead sphere:

In addition to the moment of momentum caused by the attractive force F of the respective opposing large lead spheres (distance b), the attractive force  $F_2$  of the respectively more distant sphere (distance d') gives rise to an antitorque moment (see Fig. 3). Thus, to formulate it more precisely than equation (II), the following applies for the moment of momentum  $M_1$ :

$$M_{\rm I}=2~(F+F_2^{\perp})\cdot d_{\rm I}$$

where

$$F_2^{\perp} = -F_2 \cdot \frac{b}{d'}$$

is the component of force

$$F_2 = F \cdot \frac{b^2}{{d'}^2}$$

which is acting perpendicular to the transverse beam (see (I)). The gravitational constant G is thus greater than calculated in (X) by the correction factor


$$K = \frac{F}{F + F_2^{\perp}} = \frac{1}{1 - \frac{b^3}{d'^3}}$$
(XI)

With

$$d' = \sqrt{(2d)^2 + b^2}$$

From LD Didactic Instructions P1.1.3.1 pp3

we can calculate the numerical value K = 1.083.



$$\begin{array}{ll} \mbox{Error analysis} & 1.49655 \leq M < 1.49665 \mbox{kg} \\ \hline G = \frac{\pi^2 b^2 d}{MT^2} (\theta_{\rm A} - \theta_{\rm B}) & 49.5 \leq d < 50.5 \mbox{mm} \\ & 49.5 \leq d < 50.5 \mbox{mm} \\ & 46.45 \leq b < 46.55 \mbox{mm} \\ & 46.45 \leq b < 46.55 \mbox{mm} \\ & 46.45 \leq b < 46.55 \mbox{mm} \\ & 615.58 \leq T < 616.58 \\ \hline \mbox{tan}^{-1} \bigg( \frac{1292}{6821} \bigg) - \mbox{tan}^{-1} \bigg( \frac{1046}{6819} \bigg) \leq \theta_{\rm A} - \theta_{\rm B} < \mbox{tan}^{-1} \bigg( \frac{1294}{6819} \bigg) - \mbox{tan}^{-1} \bigg( \frac{1044}{6821} \bigg) \\ \hline \mbox{G}_{\rm max} = K \times \frac{\pi^2 \big( 46.55 \times 10^{-3} \big)^2 \times 50.5 \times 10^{-3}}{1.49655 \times 615.5^2} \bigg( \mbox{tan}^{-1} \bigg( \frac{1294}{6819} \bigg) - \mbox{tan}^{-1} \bigg( \frac{1044}{6821} \bigg) \bigg) \\ \hline \mbox{We shall} \\ \mbox{assume a } +/- \mbox{precision} \\ \mbox{error of 1 mm} \\ \mbox{for all length} \\ \mbox{measurements} \\ \hline \mbox{G}_{\rm min} = K \times \frac{\pi^2 \big( 46.45 \times 10^{-3} \big)^2 \times 49.5 \times 10^{-3}}{1.49665 \times 616.5^2} \bigg( \mbox{tan}^{-1} \bigg( \frac{1292}{6821} \bigg) - \mbox{tan}^{-1} \bigg( \frac{1046}{6819} \bigg) \bigg) \\ \hline \mbox{We shall} \\ \mbox{assume a } +/- \mbox{precision} \\ \mbox{error of 1 mm} \\ \mbox{for all length} \\ \mbox{measurements} \\ \hline \mbox{G}_{\rm min} = 7.0 \times 10^{-11} \mbox{m}^3 \mbox{kg}^{-1} \mbox{s}^2 \\ \hline \mbox{So we have about a 5\% systematic offset} \\ \hline \end{tabular}$$

So final result is:

and a similar random error  $7.0 \times 10^{-11} \text{ m}^3 \text{kg}^{-1} \text{s}^{-2} < G < 7.4 \times 10^{-11} \text{ m}^3 \text{kg}^{-1} \text{s}^{-2}$